

X-ray powder diffraction – a practical guide















Dinnebier 0; InternationalTables 0; Will 0; Young 0; Pecharsky 0

X-ray hitting condensed matter




X-ray hitting condensed matter




Dinnebier Pre6

Debye-Scherrer cones from a polycrystalline sample




Bragg-Brentano and Guinier diffractometer




Bragg-Brentano diffractometer with monochromator



Generation of x-ray





Generation of x-ray



Optimum voltage ~4 times characteristic energy (~30 kV for Cu anodes)

Generation of x-ray





	Target	Filter	Incident beam* $\frac{I(K\alpha)}{I(K\beta)}$	Filter thic $\frac{I(K\alpha)}{I(K\beta)}$ in trans	$=\frac{500}{1}$	$\frac{I(K\alpha) \text{ trans.}}{I(K\alpha) \text{ incident}}$
ı				mg/cm ²	in.	
	Mo Cu Co Fe Cr	Zr Ni Fe Mn V	5.4 7.5 9.4 9.0 8.5	77 18 14 12 10	0.0046 0.0008 0.0007 0.0007 0.0006	0.29 0.42 0.46 0.48 0.49

* This is the intensity ratio at the target [G.11, Vol. 3, p. 71]. This ratio outside the x-ray tube will be changed somewhat by the differential absorption of $K\alpha$ and $K\beta$ by the tube window, typically beryllium, 0.01 inch (0.25 mm) thick.

Suppression of $K\beta$ radiation by filter with lighter neighbor element in periodic table



Well prepared samples at the right sample holder is the key for success!!!



Hygiene in preparing the powder is the second key for success!!!




Fig. 9–1 Back-reflection pinhole patterns of recrystallized aluminum specimens; grain size decreases in the order (a), (b), (c), (d). Filtered copper radiation.




Fig. 9-3 Changes in hardness and diffraction lines of 70 Cu-30 Zn specimens, reduced in thickness by 90 percent by cold rolling, and annealed for 1 hour at the temperatures indicated in (a). (b), (c), and (d) are portions of back-reflection pinhole patterns of specimens annealed at the temperatures stated (filtered copper radiation).


Bragg-Brentano diffractometer for the desk




Goniometer & optics (incident)



Goniometer & optics (receiving)



Example: growth of PrAuSi out of Sn flux



Phase analysis with the PDF database

I-0024 MAJOR CORRECTION

POWDER DIFFRACTION FILE

Sets 1-5 (Revised)

Inorganic Volume, No. PD1S-5iRB

Published by the

JOINT COMMITTEE ON POWDER DIFFRACTION STANDARDS
1601 Park Lane, Swarthmore, Pennsylvania 19081
ILS A

33-1161

0068 d 1-0024	12.8	3.60	6.9	12.8	ZeOC	£2.8H20					
1/I ₁ 1-0024	100	83	67	100		ONIUM O		DE			
Rad, W	hKa	A 0.709		Filter ZnOe		ďΑ	I/I,	hki	4 Y	I/I ₁	hkl
	6 INCHES			Coll.		12.8	100	100	2.07	13	
	AL IBRATEI			corr. abs.?	No	10.6	27	001	2.00	7	
Ref. H		DUNIFO				7.9	20	101	1.91	13	
						6.9	67	111	1.81	13	
C T	ETRAGONAL		S.C.			4.80	13	102	1.71	13	ì
Ba 12.	9 b.		LO.6 A	c				1 1			
a 10.	, 04		2	·		4.12	27	221	1.62	13	
	.P. Cons	UTED FROM	a powhrp t	1474		3.82	20	212,311		7	1
IUI.				,		3.60	83	320	1.51	7	1
						3.24	40	400	1.46	7 7	ı
ŧα			552f y 1.5	53 Sign	+	2.96	7	213	1.42	7	1
2V	D	mp	Color			2.74	7	402		l	1
Ref. W						2.55	7	430, 104	1	}	ì
						2.39	1 7	520		1	1
						2.39	13	304,441	ı	1	1
						2.15	20	600	l .		1
						4.15	۳.	INDEXED		1	1
						1		BY B.P.			
									1		

,,	-1101										33-1162
d	3.34	4.26	1.82	4.26	SiO ₂						*
1/11	100	22	14	22	Silicon Oxide				Quartz	low.	
0.4	C.V. 1	1 540500	Filter M	ono, Dia		d A	1/11	hkt	d A	· 1/1 ₁	hk1
Cut of	ff	1/1, Di	ffractome	ter 1/1		4.257 3.342 2.457 2.282	22 100 8 8	100 101 110 102	1.2285 1.1999 1.1978 1.1843	1 2 1 3	220 213 221 114
Sys. 1 a ₀ 4.	Hexagonal 9133(2)b _ξ β Ibid.			P3 ₁ 21 3(4) A Z 3		2.237 2.127 1.9792 1.8179	6 4	111 200 201 112	1.1804 1.1532 1.1405 1.1143	3 1 <1 <1 <1	310 311 204 303
ea 2V Ref.		nωβ 1. 2.656	544 ¢		Sign + rColorless	1.8021 1.6719 1.6591 1.6082	<1 4 2 <1	003 202 103 210	1.0813 1.0635 1.0476 1.0438	1	312 400 105 401
Sampl of St	e from th	ground si	ngle crys	tals of o	onal Bureau ptical quality,	1.5418 1.4536 1.4189	1 <1	211 113 300	1.0347 1.0150 0.9898	1	214 223 402
Silic Fin =	on (ac=S.	43088Å) u 126,31).	ern at 25 sed as in Quartz g	ternal st	andard.	1.3820 1.3752 1.3718 1.2880	8 2	212 203 301 104	.9873 .9783 .9762 .9636	1 <1	313 304 320 205
						1.2558	2	302	6 refle	tions	to 0.9089

0831 d 1-0363	4.30	3.81	4.08	4.30	Sic) ₂			S102			
I/I ₁ 1-0378	100	67	33	100	Sii	LICON DI	Oxide		TAIDYM	117E (LOI	FORM)	
Rad N	hka	1 0,709		Filter ZnOo		фÅ	I/I ₁	hkl	d Å	I/I ₁	PRI	
	6 INCHES			Coll.		4.30	100		1.44	3		
	ALIBRATES		d	corr. abs.?	No	4.08	33		1.40	7		
Ref. H	1					3.61	67		1.36	3		
						3.43	1		1.31	4		
See C	RTHORHOME	110	S.G.			3.21	1		1.25	3	!	
		7.1 ce		.578 C C	954	ı					1	
a a	В	Ý		64 ?		2.96	17		1.19	5		
Ref. W	r (FOR OR	ILY ONE TY	PE OF TRI	DYMITE)		2.80	3 27		1.15	3	1	
						2.49	11		1.10	1		
8 a 1	.478	D 0 B 1.4	79 E v 1.4	81 Slen		2.08	11 5					
ZV 35	D 2.2	26 mp	Color C	OLORLESS		2.00	,					
Ref. W	N, C.C. (FOR ONLY	ONE TYPE	OF TRIDYM	ITE)	1.84	3					
						1.69	8					
T.P.	το β1 (LC	WER HIGH-	TRIDYUITE) AT 117	C	1.64	4					
T.P.	το β ₂ (υρ	PER HIGH-	TRIDYMITE	C	1.60	4	1		1	ı		
		AND REFR				1.53	5					
		OCCURS II				1					ı	
		FLEOTION				/				<u>'</u> .\ e	· F	

Phase analysis with the PDF database

EXPLANATION OF THE FORMAT

· 1	a 1 b	1 c			7				8	
100		3			-	IA.		**	u.	-
i i		- 5 ^{''}	C	Sep		-	9			
		6								

The card shown above has spaces numbered from 1 to 10 inclusive. The explanation of the symbols in the various spaces is as follows:

Spaces la 1h and 1c

These contain the interplanar spacings corresponding to the three strongest lines in the diffraction pattern, chosen from the forward reflections, where 24<90° Note that where $2\theta = 90^{\circ}$, the value of dA depends upon the wavelength of the radiation used

Space 1d

This gives the largest interplanar spacing found for the specimen.

Spaces 2a, 2b, 2c and 2d

These contain the intensities of the lines in Spaces 1a, 1b, 1c, 1d, referred to the strongest line as 100. This intensity occasionally is given a number greater than 100 when it is very much stronger than the rest of the pattern.

Space 3

Rad.-Source of the x-rays (Mo, Cu, etc.)

λ-The wave length of the x-rays used in Angstroms.

Filter—The substance used to filter out extraneous wave lengths.

Dia.-The diameter of the cylindrical film holder.

Cut off-The longest spacing possible with the apparatus used.

I/I.-The method used to measure relative intensities. This ordinarily means either diffractometer, calibrated strips of photographic film or visual inspection. The relative intensities are expressed on a scale 0-100 in steps of 10 for visually estimated values, and more exactly for measured values.

Ref.-Source of the data listed in Spaces 3 and 9.

Space 4

Sys.-Crystallographic system to which the specimen belongs.

S. G.-Space group, listed according to the Schoenflies system and the system given in "International Tables for X-ray Crystallography" (1952).

ao, bo, and co-Lattice parameters. A = ao/bo C = co/bo.

a. R. v-Interaxial angles.

Z-The numbers of chemical formula units per unit of structure. For chemical elements. Z represents the number of atoms per unit of structure; for compounds. Z represents the number of formula units per unit cell. Spaces 7 and 8 show usually the simplest atomic formula for the compound. The "dot" formula has been multiplied by an appropriate factor to make it agree with the formula in Space 7.

Dx-Density calculated from x-ray measurements Ref.-Source of the data listed in Space 4

Space 5

ea, nωβ, and ey-Indices of refraction.

Sign and 2V have the customary crystallographic meanings.

D-Measured density

mn-Melting point

Color-The color of the specimen as ordinarily seen or as seen when examined by microscopic methods. Occasionally other data are listed in this space, such as hardness (H) and luster of minerals. Ref.-Source of the data listed in Space 5.

Space 6

This contains further pertinent information, such as chemical analysis of the specimen, source of the sample, heat treatment, temperature at which the nattern was made etc.

Space 7

Chemical formula and name of the specimen. The formula may be omitted in cases of too complex compositions.

Space 8

"Dot" or structural formula for the specimen, when available, above the mineralogical or common name, if any, of the specimen. Parentheses around the name indicates a synthetic material. A * in the upper right corner of this space indicates that the card contains data of high reliability: a O indicates low reliability.

Space 9

This provides columns of interplanar spaces, relative intensities, and Miller indices.

The following abbreviations may be used in Space 9:

b = Broad, fuzzy or diffuse line

d - Doublet

n = Line not given by all sources

nc = Line not accounted for by the proposed unit cell

ni = Line cannot be indexed with given unit cell

np = Index not permitted by given space group

R = Intensity uncertain owing to presence of, or overlapping of, R lines tr = Trace

+ = Additional indices are possible

Space 10

The location of the identification number of the card. Because of changes in the d-spacings of many of the patterns in these sets, the identification number no longer shows the Hanawalt order of the card in the File and should be considered an arbitrary number only. PDF DataBase 5: 6

Phase analysis with the PDF database

24.7x

15.1, 12.2, 11.3_x

4.94, 3.54, 3.10, 2.48,

12.26 8.65, 7.85,

12.1, 10.6, 5.79, 4.35,

7.44

9.80, 8.57,

11.1

4.25.

Powder Diffraction File Search Manual



Published by the

INTERNATIONAL CENTRE FOR DIFFRACTION DATA

1601 PARK LANE • SWARTHMORE, PA 19081-2389 • U.S.A

FINK

POWDER	DIFFRAC	TION	FILE
RET	RIEVAL II	NDEX	

for

INORGANIC COMPOUNDS 1972

Publication PDIS-22f

Published by the

JOINT COMMITTEE ON POWDER DIFFRACTION STANDARDS 1601 Park Lane, Swarthmore, Pennsylvania 19081

•	30.06	15.0 _e	9.701	4.53x 7.30g	3.293	2.983	8.57 2.54 3.27 2.98	1.51,	Hf (OH) 2504.2H20 Na-K-Ca-Mg-Al-Si-O-OH-H2O Mg-kl3514020 (OH) 10-4820 Ca2Cu9 (ASO4) 4 (OH) 10.10H2O	12- 231 13- 190	1-134-C 3 1-46-8 5 1-52-E 4 1-80-F 7
	23.1 _x 16.1 ₄ 14.5 _x	13.9 ₁ 13.8 ₄ 12.9 ₁	11.9 _x 10.6 _s 11.8 ₁	9.70 ₁ 9.07 ₅ 7.20 ₁	9.10 ₁ 7.12 _x 4.88 ₁	8.00 ₁ 4.82 ₇ 4.25 ₁	2.83 ₃ 6.91 ₁ 4.72 ₆ 2.60 ₁ 2.56 ₅	3.18 ₁ 3.83 ₂ 0.00 ₁	Hg Fe Alsi Pe ₄ (PO ₄) ₃ (OB) ₃ . 12H ₂ O C ₂ H ₈ N ₂ O ₆ U.3 H ₂ O HaGa ₄ O ₆ Al-7i-Fe-K-Na-Si-OH	14- 331 19- 69 21-1114	I-109-F 1 I-59-8 6 I-103-E 2 I-141-B 3 I- 26-C 3
	05.5		2 20							12- 222	

3.33 2.37

5.98

2.90.

10.0 and over

Al-Pe-Hg-Ca-Na-K-Si-O.H.O

Na-Ca-K-H-O-A1-Si-O-R.38+0

(Al, Fe) , AsO, (OH) , .5H, O

Rf (OH) 2 SO4 . 2H2O

V205.3820


Pile No. Piche No.

14- 183 I- 57-E12


11- 146 T- 39-811

21- 361 I-134-C 3 7- 332 I- 26-C 4


Example: growth of PrAuSi out of Sn flux




If your phase is not in the database – search for isostructural compounds...



The best tool to calculate diffraction pattern, to verify structure data and more... $_{\tiny Powdercell\ 1}$



The best tool to calculate diffraction pattern, to verify structure data and more...



The best tool to calculate diffraction pattern, to verify structure data and more...

Powdercell 3

① Cmm2

 $C_{2\nu}^{11}$

mm2


Orthorhombic

2 No. 35

Cmm2

Patterson symmetry Cmmm





(3)





- (4) Origin on mm2
- (5) Asymmetric unit $0 \le x \le \frac{1}{4}$; $0 \le y \le \frac{1}{2}$; $0 \le z \le 1$
- (6) Symmetry operations

For (0,0,0)+ set (1) 1

(1) 1 (2) 2 0,0,z (3)

(3) m x,0,z (4) m 0,y,z

For $(\frac{1}{2}, \frac{1}{2}, 0)$ + set (1) $t(\frac{1}{2}, \frac{1}{2}, 0)$

(2) 2 1,1,z

(3) $a = x, \frac{1}{2}, z$ (4) $b = \frac{1}{2}, y, z$

Headline: Section 2.2.3.
 Short Hermann-Mauguin symbol Schoenflies symbol Crystal class (Point group) Crystal system (Section 2.2.4 and Chapter 12.2) (Chapters 12.1 and 12.2) (Section 10.1.1 and Chapter 12.1) (Section 2.1.2)

Number of space group Full Hermann–Mauguin symbol Patterson symmetry
[Same as in IT (1952)] (Section 2.2.4 and Chapter 12.3) (Section 2.2.5)

- 3. Space-group diagrams, consisting of one or several projections of the symmetry elements and one illustration of a set of equivalent points in general position. The numbers and types of the diagrams depend on the crystal system. The diagrams and their axes are described in Section 2.26 the graphical symbols of symmetry elements are listed in Chapter 1.4. For monoculinit source remuse see Section 2.2.16 for orthorhombics estimes see Section 2.2.6.4.
- (4) Origin of the unit cell: Section 2.2.7. The site symmetry of the origin and its location with respect to the symmetry elements are given.
- (5) Asymmetric unit: Section 2.2.8. One choice of asymmetric unit is given.
- (6) Symmetry operations: Section 2.2.9 and Part 11. For each point \(\tilde{x}, \tilde{y}, \tilde{z}\) of the general position that symmetry operation is listed which transforms the initial point x, y, \(\tilde{z}\) into the point under consideration. The symbol describes the nature of the operation, its glide or screw component (given between parentheses), if present, and the location of the corresponding symmetry element.

The symmetry operations are numbered in the same way as the corresponding coordinate triplets of the general position. For centred space groups the same numbering is applied in each block, $e_i e_i$ under 'For $(k_i, k_i, 0)$ + set'.

[Continued on inside back cover]

1.4.1. Symmetry planes normal to the plane of projection (three dimensions) and symmetry lines in the plane of the figure (two dimensions)

Symmetry plane or symmetry line	Graphical symbol	Glide vector in units of lattice translation vectors parallel and normal to the projection plane	Printed symbol
Reflection plane, mirror plane Reflection line, mirror line (two dimensions)		None	m
'Axial' glide plane Glide line (two dimensions)		$\frac{1}{2}$ lattice vector along line in projection plane $\frac{1}{2}$ lattice vector along line in plane	a, b or c
'Axial' glide plane		$\frac{1}{2}$ lattice vector normal to projection plane	a, b or c
'Double' glide plane* (in centred cells only)		Two glide vectors: $\frac{1}{2}$ along line parallel to projection plane, $\frac{1}{2}$ normal to projection plane	e
'Diagonal' glide plane		One glide vector with two components: \[\frac{1}{2} \] along line parallel to projection plane, \[\frac{1}{2} \] normal to projection plane	n
'Diamond' glide plane† (pair of planes; in centred cells only)		\$\frac{1}{4}\$ along line parallel to projection plane, combined with \$\frac{1}{4}\$ normal to projection plane (arrow indicates direction parallel to the projection plane for which the normal component is positive)	d

^{*} For further explanations of the 'double' glide plane e see Note (iv) below and Note (x) in Chapter 1.3.

[†] See footnote § to Section 1.3.1.

1.4.2. Symmetry planes parallel to the plane of projection

Symmetry plane	Graphical symbol*	Glide vector in units of lattice translation vectors parallel to the projection plane	Printed symbol
Reflection plane, mirror plane		None	m
'Axial' glide plane		$\frac{1}{2}$ lattice vector in the direction of the arrow	a, b or c
'Double' glide plane† (in centred cells only)	\	Two glide vectors: $\frac{1}{2}$ in either of the directions of the two arrows	e
'Diagonal' glide plane		One glide vector with <i>two</i> components $\frac{1}{2}$ in the direction of the arrow	n
'Diamond' glide plane‡ (pair of planes; in centred cells only)	3 8	$\frac{1}{2}$ in the direction of the arrow; the glide vector is always half of a centring vector, <i>i.e.</i> one quarter of a diagonal of the conventional face-centred cell	d

^{*}The symbols are given at the upper left corner of the space-group diagrams. A fraction h attached to a symbol indicates two symmetry planes with 'heights' h and $h + \frac{1}{2}$ above the plane of projection; e.g. $\frac{1}{8}$ stands for $h = \frac{1}{8}$ and $\frac{1}{8}$. No fraction means h = 0 and $\frac{1}{2}$ (cf. Section 2.2.6).

[†] For further explanations of the 'double' glide plane e see Note (iv) below and Note (x) in Chapter 1.3.

[‡] See footnote § to Section 1.3.1.

1.4.5. Symmetry axes normal to the plane of projection and symmetry points in the plane of the figure

Symmetry axis or symmetry point	Graphical symbol*	Screw vector of a right-handed screw rotation in units of the shortest lattice translation vector parallel to the axis	Printed symbol (partial elements in parentheses)
Identity	None	None	1
Twofold rotation axis Twofold rotation point (two dimensions)	•	None	2
Twofold screw axis: '2 sub 1'	9	1/2	21
Threefold rotation axis Threefold rotation point (two dimensions)	A	None	3
Threefold screw axis: '3 sub 1'	A	1	31
Threefold screw axis: '3 sub 2'	A	3	32
Fourfold rotation axis Fourfold rotation point (two dimensions)	+ =	None	4 (2)
Fourfold screw axis: '4 sub 1'	→ 1	1	41 (21)
Fourfold screw axis: '4 sub 2'	- - →	1	42 (2)
Fourfold screw axis: '4 sub 3'	- ∳ ≒	3 4	43 (21)
Sixfold rotation axis Sixfold rotation point (two dimensions)	•	None	6 (3,2)
Sixfold screw axis: '6 sub 1'	₩.	ļ	6, (3, 2,)
Sixfold screw axis: '6 sub 2'	•	1	62 (32, 2)
Sixfold screw axis: '6 sub 3'		1/2	63 (3,21)
Sixfold screw axis: '6 sub 4'	•	3	64 (31, 2)
Sixfold screw axis: '6 sub 5'	₩	5	65 (32, 21)
Centre of symmetry, inversion centre: '1 bar' Reflection point, mirror point (one dimension)	۰	None	ī
Inversion axis: "3 bar"	Δ	None	3 (3, 1)
Inversion axis: '4 bar'	◆ ≥	None	ā (2)
Inversion axis: '6 bar'		None	$\bar{6} \equiv 3/m$
Twofold rotation axis with centre of symmetry		None	2/m (Ī)
Twofold screw axis with centre of symmetry	9	1/2	2 ₁ /m (Ī)
Fourfold rotation axis with centre of symmetry	• •	None	4/m (4, 2, 1)
'4 sub 2' screw axis with centre of symmetry	ý <u>-</u>	1/2	$4_2/m \ (\bar{4},2,\bar{1})$
Sixfold rotation axis with centre of symmetry	•	None	6/m (6,3,3,2,1)
'6 sub 3' screw axis with centre of symmetry	6	1/2	63/m (6,3,3,21,1)

^{*} Notes on the 'heights' h of symmetry points $\bar{1}, \bar{3}, \bar{4}$ and $\bar{6}$:

⁽¹⁾ Centres of symmetry I and 3, as well as inversion points 4 and 6 on 4 and 6 axes parallel to [001], occur in pairs at 'heights' h and h + ½. In the space-group diagrams, only one fraction h is given, e.g. ½ stands for h = ½ and ½. No fraction means h = 0 and ½. In orbit space groups, however, both fractions are given for vertical 4 axes, including h = 0 and ½.

incoming n = 0 and f / m and 6/m contain vertical 4 and 6 axes; their 4 and 6 inversion points coincide with the centres of symmetry. This is not indicated in the space-group diagrams.

⁽³⁾ Symmetries 4:/m and 6:/m also contain vertical 4 and 6 axes, but their 4 and 6 inversion points alternate with the centres of symmetry; Le. I points at h and h + ½ interleave with 4 or 6 points at h + 2 and h + ½ interleave with 4 or 6 points at h + 2 and h + ½. In the tetragonal and becaugeousl spoce-group diagrams, only one fraction for 1 and one for 4 or 6 is given. In the cubic diagrams, after fractions are listed for 4 n/n; e.g. a Price Noc. 2023; 1: 0. (e. 4: ½.).


1.4.6. Symmetry axes parallel to the plane of projection

Symmetry axis	Graphical symbol	*	Screw vector of a right-handed scr rotation in units of the shortest latt translation vector parallel to the ax	ice (partial elements
Twofold rotation axis	+ +	† †	None	2
Twofold screw axis: '2 sub 1'	~ ~	1 4	$\frac{1}{2}$	21
Fourfold rotation axis	1-1	T	None	4 (2)
Fourfold screw axis: '4 sub 1'	1 4	7 🖈	To sdi	41 (21)
Fourfold screw axis: '4 sub 2'	JJ	→	sdnoza -12	4 ₂ (2)
Fourfold screw axis: '4 sub 3'	⊮ - ∦	+ 4) 2 2 3 4 3 4 4 4 4 4 4 4 4 4 4 4 4 4 4 4	4 ₃ (2 ₁)
Inversion axis: '4 bar'	88	7 1	none None	<u>4</u> (2)
Inversion point on '4 bar'-axis	₽	4	ii -	4 point

^{*} The symbols for horizontal symmetry axes are given outside the unit cell of the space-group diagrams. Twofold axes always occur in pairs, at 'heights' h and $h + \frac{1}{2}$ above the plane of projection; here, a fraction h attached to such a symbol indicates two axes with heights h and $h + \frac{1}{2}$. No fraction stands for h = 0 and $\frac{1}{2}$. The rule of pairwise occurrence is not valid for the horizontal fourfold axes in cubic space groups; here, all heights are given, including h = 0 and $\frac{1}{2}$. This applies also to the horizontal $\frac{1}{4}$ axes and the $\frac{1}{4}$ inversion points located on these axes.

System	Axial lengths and angles	Bravais lattice	Lattice symbol
Cubic	Three equal axes at right angles $a=b=c$, $\alpha=\beta=\gamma=90^{\circ}$	Simple Body-centered Face-centered	P I F
Tetragonal	Three axes at right angles, two equal $a=b\neq c$, $\alpha=\beta=\gamma=90^{\circ}$	Simple Body-centered	P
Orthorhombic	Three unequal axes at right angles $a \neq b \neq c$, $\alpha = \beta = \gamma = 90^{\circ}$	Simple Body-centered Base-centered Face-centered	P I C F
Rhombohedral*	Three equal axes, equally inclined $a=b=c$, $\alpha=\beta=\gamma\neq90^{\circ}$	Simple	R
Hexagonal	Two equal coplanar axes at 120°, third axis at right angles $a=b\neq c$, $\alpha=\beta=90^{\circ}$, $\gamma=120^{\circ}$	Simple	Р
Monoclinic	Three unequal axes, one pair not at right angles $a \neq b \neq c$, $\alpha = \gamma = 90^{\circ} \neq \beta$	Simple Base-centered	P C
Triclinic	Three unequal axes, unequally inclined and none at right angles $a \neq b \neq c$, $\alpha \neq \beta \neq \gamma \neq 90^{\circ}$	Simple	Р

^{*} Also called trigonal.





ORTHORHOMBIC ORTHORHOMBIC



(F)



(P)

(R)

① CONTINUED No. 35 Cmm2

Reflection conditions General: $\begin{array}{ll} \text{Bis} : h + k = 2n \\ \text{bis} : h + k = 2n \\ \text{bis} : h = 2n \\ \text{fiol} : h = 2n \\ \text{Special: as above, plus no extra conditions} \\ \text{no extra co$

(2) Generators selected (1): t(1,0,0); t(0,1,0); t(0,0,1); $t(\frac{1}{2},\frac{1}{2},0)$; (2); (3)

(3)	Positions

	Multiplicity,		Coordinates				
Wyckoff letter, Site symmetry				(0,0,0)+	$({\scriptstyle\frac{1}{2}},{\scriptstyle\frac{1}{2}},0)+$		
8	f	1	$(1)\;x,y,z$	(2) \bar{x}, \bar{y}, z	$(3)\ x,\bar{y},z$	(4) \bar{x}, y, z	
4	e	m	0, y, z	$0, \bar{y}, z$			
4	d	$.\ m\ .$	x, 0, z	$\vec{x}, 0, z$			
4	С	2	$\frac{1}{4}, \frac{1}{4}, z$	$\frac{1}{4}$, $\frac{3}{4}$, z			
2	Ь	m m 2	$0, \frac{1}{2}, z$				
2			0.0 -				

(4) Symmetry of special projections

Along [001] c 2 mm	Along [100] p1m1	Along [010] p 1 1 m
$\mathbf{a}' = \mathbf{a}$ $\mathbf{b}' = \mathbf{b}$	$\mathbf{a}' = \frac{1}{2}\mathbf{b}$ $\mathbf{b}' = \mathbf{c}$	$\mathbf{a}' = \mathbf{c}$ $\mathbf{b}' = \frac{1}{2}\mathbf{a}$
Origin at 0, 0, z	Origin at x, 0, 0	Origin at 0, y, 0

(5) Maximal non-isomorphic subgroups

	[2] Cm 11 (Cm, 8)	(1; 4)+
	[2] C112 (P2, 3)	(1; 2)+
IIa	[2] Pba2 (32)	1; 2; $(3; 4) + (\frac{1}{2}, \frac{1}{2}, 0)$
	[2] Pbm2 (Pma2, 28)	1; 3; (2; 4) + (1, 1, 0)
	[2] Pma2 (28)	1; 4; $(2; 3) + (\frac{1}{2}, \frac{1}{2}, 0)$
	[2] Pmm2 (25)	1; 2; 3; 4

IIb [2] Ima2 (c' = 2c) (46); [2] Ibm2 (c' = 2c) (Ima2, 46); [2] Iba2 (c' = 2c) (45); [2] Imm2 (c' = 2c) (44); [2] Ccc2 (c' = 2c) (37); [2] Cmc2, (c' = 2c) (36); [2] Ccm2, (c' = 2c) (Cmc2, 36)

(6) Maximal isomorphic subgroups of lowest index

He [2] Cmm2 (c' = 2c) (35); [3] Cmm2 (a' = 3a or b' = 3b) (35)

(7) Minimal non-isomorphic supergroups

- I [2] Cmmm (65); [2] Cmme (67); [2] P4mm (99); [2] P4bm (100); [2] P4,cm (101); [2] P4,nm (102); [2] $P\bar{4}2m$ (111); [2] $P\bar{4}2,m$ (113); [3] P6mm (183)
- II [2] Fmm2 (42); [2] Pmm2 ($\mathbf{a}' = \frac{1}{2}\mathbf{a}, \mathbf{b}' = \frac{1}{2}\mathbf{b}$) (25)

- Headline in abbreviated form.
- ② Generators selected: Sections 2.2.10 and 8.3.5. A set of generators, as selected for these Tables, is listed in the form of translations and numbers of general-position coordinates. The generators determine the sequence of the coordinate triplets in the general position and of the corresponding symmetry operations.
- 3 Positions: Sections 2.2.11 and 8.3.2. The general Wyckoff position is given at the top, followed downwards by the various special Wyckoff positions with decreasing multiplicity and increasing site symmetry. For each general and special position is multiplicity, Wyckoff letter, oriented site-symmetry symbol, as well as the appropriate coordinate triplets and the reflection conditions, are listed. The coordinate triplets of the general position are numbered sequentially: cf. Symmetry operations.

Oriented site-symmetry symbol (third column): Section 2.2.12. The site symmetry at the points of a special position is given in oriented form.

Reflection conditions (right-most column): Section 2.2.13.

[Lattice complexes are described in Part 14; Tables 14.2.3.1 and 14.2.3.2 show the assignment of Wyckoff positions to Wyckoff sets and to lattice complexes.]

- Symmetry of special projections: Section 2.2.14. For each space group, orthographic projections along three (symmetry) directions are listed. Given are the projection direction, the plane group of the projection, as well as the axes and the origin of the projected cell.
- (5) Maximal non-isomorphic subgroups: Sections 2.2.15 and 8.3.3.

Type I: translationengleiche or t subgroups;

Type Ha: klassengleiche or k subgroups, obtained by 'decentring' the conventional cell; applies only to space groups with centred

Type IIb: klassengleiche or k subgroups, obtained by enlarging the conventional cell.

Given are:

For types I and IIa: Index [between brackets]; 'unconventional' Hermann-Mauguin symbol of the subgroup; 'conventional' Hermann-Mauguin symbol of the subgroup, if different (between parentheses); coordinate triplets retained in subgroup. For type IIb: Index [between parentheses] coordinate triplets retained in subgroup. Both packets: 'unconventional' Hermann-Mauguin symbol of the subgroup: basis-vector relations between

For type IIb: Index [between brackets]: "unconventional" Hermann-Mauguin symbol of the subgroup; basis-vector relations between group and subgroup (between parentheses); "conventional" Hermann-Mauguin symbol of the subgroup, if different (between parentheses).

(6) Maximal isomorphic subgroups of lowest index: Sections 2.2.15, 8.3.3 and 13.1.2.

Type IIc: klassengleiche or k subgroups of lowest index which are of the same type as the group, i.e. have the same standard Hermann–Mauguin symbol. Data as for subgroups of type IIb.

Minimal non-isomorphic supergroups: Sections 2.2.15 and 8.3.3.
The list contains the reverse relations of the subgroup tables; only types I (t supergroups) and II (k supergroups) are distinguished Data as for subgroups of type IIb.

Problems describing a structure - Rhombohedral unit cell




Fig. 5.1.3.6. Unit cells in the rhombohedral lattice: same origin for all cells. The basis of the hombohedral cell is labelled a, b, c. Two settings of the triple hexagonal cell are possible with respect to a primitive rhombohedral cell: The *obverse setting* with the lattice points $0, 0, 0, \frac{2}{3}, \frac{1}{3}, \frac{1}{3}, \frac{1}{3}, \frac{1}{3}$ has been used in *International Tables* since 1952. Its general reflection condition is -h + k + l = 3n. The *reverse setting* with lattice points $0, 0, 0, \frac{1}{3}, \frac{1}{3}, \frac{1}{3}$ was used in the 1935 edition. Its general reflection condition is h + k + l = 3n. (a) Obverse setting of triple hexagonal cell a, b, c, b. (c) Primitive rhombohedral cell a, b, c. (b) Reverse setting of triple hexagonal cells in relation to the primitive rhombohedral cell (---) lower edges), a, b, c in relation to the three triple hexagonal cells in obverse setting a_1, b_1, c' ; a_2, b_2, c' ; a_3, b_3, c' . Projection along c'. (d) Primitive rhombohedral cell (---) lower edges), a, b, c in relation to the three triple hexagonal cells in reverse setting a_1, b_1, c' ; a_2, b_2, c' ; a_3, b_3, c' . Projection along c'. (a) Primitive rhombohedral cell (---) lower edges), a_3, b_3, c' in relation to the three triple hexagonal cells in reverse setting a_1, b_1, c' ; a_3, b_3, c' . Projection along c'.

Problems describing a structure - Rhombohedral unit cell

R3

3

R3

RHOMBOHEDRAL AXES


No. 146

R3

No. 146

R3

HEXAGONAL AXES



Generators selected (1); t(1,0,0); t(0,1,0); t(0,0,1); $t(\frac{1}{2},\frac{1}{2},\frac{1}{2})$; (2)

Positions Site symmetry

Multiplicity, Coordinates Wyckoff letter.

(2) $\bar{y}, x - y, z$

 $(0,0,0)+(\frac{2}{3},\frac{1}{3},\frac{1}{3})+(\frac{1}{3},\frac{2}{3},\frac{2}{3})+$

3 a 3. 0.0.z

Symmetry of special projections

(1) x, y, z

Along [001] p3 $a' = \frac{1}{3}(2a + b)$ $\mathbf{b}' = \frac{1}{3}(-\mathbf{a} + \mathbf{b})$ Origin at 0.0.z

Along [100] p1 $\mathbf{a}' = \frac{1}{2}(\mathbf{a} + 2\mathbf{b})$ Origin at x, 0, 0

(3) $\bar{x} + y, \bar{x}, z$

 $\mathbf{b}' = \frac{1}{2}(-\mathbf{a} - 2\mathbf{b} + \mathbf{c})$

Generators selected (1): t(1,0,0): t(0,1,0): t(0,0,1): (2)

Positions

3 b

Multiplicity. Wyckoff letter Site symmetry

(2) z, x, y

(3) y, z, x

Coordinates

X, X, X

Symmetry of special projections

Along [111] p3 $\mathbf{a}' = \frac{1}{2}(2\mathbf{a} - \mathbf{b} - \mathbf{c})$ Origin at x, x, x

 $\mathbf{b}' = \frac{1}{2}(-\mathbf{a} + 2\mathbf{b} - \mathbf{c})$

(1) x, y, z

Along [110] p1 $\mathbf{a}' = \frac{1}{2}(\mathbf{a} + \mathbf{b} - 2\mathbf{c})$ Origin at $x, \bar{x}, 0$

 $\mathbf{b}' = \mathbf{c}$

Problems describing a structure – Origin of cell

$$D_{4h}^7$$

P4/nmm


No. 129

 $P 4/n 2_1/m 2/m$

No. 129

 $P 4/n 2_1/m 2/m$

ORIGIN CHOICE 1




Origin at 4m2 at 4/nm2/g, at $-\frac{1}{4}, \frac{1}{4}, 0$ from centre (2/m)

Asymmetric unit $0 \le x \le \frac{1}{2}$; $0 \le y \le \frac{1}{2}$; $0 \le z \le \frac{1}{2}$; $y \le \frac{1}{2} - x$

Symmetry operations

- (1) 1(5) $2(0,\frac{1}{2},0)$ $\frac{1}{4},y,0$ (9) 1 1.1.0 (13) m x, 0, z
- (2) 2 0,0,z(6) $2(\frac{1}{2},0,0)$ $x,\frac{1}{4},0$
- (10) $n(\frac{1}{5}, \frac{1}{5}, 0) x, y, 0$ (14) m = 0, v, z
- (3) 4^+ $0, \frac{1}{3}, z$ (15) $m x + \frac{1}{3}, \bar{x}, z$
- (7) 2 x, x, 0(11) $\bar{4}^+$ 0.0,z; 0.0,0

ORIGIN CHOICE 2



Origin at centre (2/m) at n2, (2/m, 2, /g), at $\frac{1}{2}$, $-\frac{1}{2}$, 0 from 4m2

Asymmetric unit $-\frac{1}{4} \le x \le \frac{1}{4}; \quad -\frac{1}{4} \le y \le \frac{1}{4}; \quad 0 \le z \le \frac{1}{5}; \quad x \le y$

Symmetry operations

(13) $m = x, \frac{1}{4}, z$

- (1) 1(2) $2 + \frac{1}{2}, \frac{1}{2}, z$ (5) $2(0,\frac{1}{2},0)$ 0,y,0(9) Ī 0.0.0
 - (6) $2(\frac{1}{2},0,0)$ x,0,0(10) $n(\frac{1}{3}, \frac{1}{3}, 0)$ x, y, 0(14) $m^{-\frac{1}{4}}, y, z$
- (3) 4⁺ ½, ½, z (7) $2(\frac{1}{2}, \frac{1}{2}, 0)$ x, x, 0
- (11) $\bar{4}^{+}$ $\frac{1}{4}$, $-\frac{1}{4}$, z; $\frac{1}{4}$, $-\frac{1}{4}$, 0 (15) $m x + \frac{1}{2}, \bar{x}, z$

Problems describing a structure – Origin of cell

P4/nmm

 D_{4h}^7

P4/nmm

No. 129

 $P 4/n 2_1/m 2/m$

No. 129

 $P 4/n 2_1/m 2/m$

ORIGIN CHOICE 1

. . m

 $\bar{4} m 2$

(1) x, y, z

 $x, x + \frac{1}{2}, z$

 $\bar{x} + \frac{1}{3}, x, \bar{z}$

Positions Multiplicity Wyckoff letter Site symmetry

16 k 1

Coordinates

(3) $\bar{v} + \frac{1}{2} \cdot x + \frac{1}{2} \cdot z$ (2) \bar{x}, \bar{y}, z (6) $x + \frac{1}{3}, \vec{v} + \frac{1}{3}, \vec{z}$ (7) v, x, \bar{z}

 $\bar{x}.x + \frac{1}{2}.7$

 $x + \frac{1}{3}, x, \bar{z}$

 $x, \bar{x} + \frac{1}{2}, 7$

 $\bar{x} + \frac{1}{3}, \bar{x}, \bar{z}$

(5) $\bar{x} + \frac{1}{2}, v + \frac{1}{2}, \bar{z}$ (9) $\bar{x} + \frac{1}{3}, \bar{y} + \frac{1}{3}, \bar{z}$ (10) $x + \frac{1}{3}, y + \frac{1}{3}, \overline{z}$ (11) v, \bar{x}, \bar{z} (13) x, \bar{y}, z (14) \bar{x}, y, z (15) $\bar{y} + \frac{1}{3}, \bar{x} + \frac{1}{3}, z$

 $\bar{x}.\bar{x} + \frac{1}{2}.z$

 $x + \frac{1}{2}, \bar{x}, \bar{z}$

ORIGIN CHOICE 2

Positions Multiplicity. Wyckoff letter. Site symmetry k 1

j ..m

8 g ..2

4 m m

 $\bar{4}$ m 2

 $\bar{4} m 2$

16

(1)
$$x, y, z$$

(5) $\bar{x}, y + \frac{1}{2}, \bar{z}$
(9) $\bar{x}, \bar{y}, \bar{z}$
(13) $x, \bar{y} + \frac{1}{2}, z$

x, x, z

 $\frac{1}{4}$, y, z

 $\frac{3}{7}, \frac{1}{7}, 0$

 $\bar{x}, x + \frac{1}{2}, \bar{z}$

(2)
$$\bar{x} + \frac{1}{2}, \bar{y} + \frac{1}{2}, z$$

(6) $x + \frac{1}{2}, \bar{y}, \bar{z}$
(10) $x + \frac{1}{2}, y + \frac{1}{2}, \bar{z}$
(14) $\bar{x} + \frac{1}{2}, y, z$

 $\bar{x} + \frac{1}{3}, \bar{x} + \frac{1}{3}, z$

 $x + \frac{1}{3}, \bar{x}, \bar{z}$

 $\frac{1}{7}$, $\bar{y} + \frac{1}{7}$, z

Coordinates

(3)
$$\bar{y} + \frac{1}{2}, x, z$$

(7) $y + \frac{1}{2}, x + \frac{1}{2}, \bar{z}$
(11) $y + \frac{1}{2}, \bar{x}, \bar{z}$
(15) $\bar{y} + \frac{1}{2}, \bar{x} + \frac{1}{2}, z$

 $y, \frac{1}{4}, z$

 $\bar{x}, \bar{x} + \frac{1}{2}, \frac{1}{2}$

 $x, x + \frac{1}{3}, \frac{1}{3}$

 $\bar{x}.\bar{x} + \frac{1}{2}.0$

 $x, x + \frac{1}{2}, 0$

 $x, \bar{x} + \frac{1}{2}, z$

 \bar{X} , \bar{X} , \bar{Z}

 $\bar{x} + \frac{1}{2}, x, z$

 $\bar{v} + \frac{1}{2}, \frac{1}{2}, z$


 $x + \frac{1}{2}, x + \frac{1}{2}, \bar{z}$

 $\frac{1}{2}, \frac{1}{2}, 0$

0.0.0

 $\frac{3}{7}, \frac{3}{7}, \overline{Z}$

 $\frac{1}{4}, \frac{3}{4}, 0$



Extract the reflection list: (hkl) – position - intensity

Preparation of a reflection list for further use

- 1. Collect a diffraction pattern from the pure phase. If pure phase material is not available, then the phase should constitute the bulk of the sample.
- 2. Run the Rietveld program in the LeBail fitting mode using the assigned space group and unit cell parameters. From the refined list of intensities, create a file containing h, k, l, M, d, 2θ and I, where h, k and l are the Miller indices of the reflection, M is the reflection multiplicity, d is the d-spacing of the reflection, 2θ is the Bragg angle and I is the reflection intensity.
- 3. Depending on which Rietveld program has been used, it might be necessary to remove the effect of the Lorentz-polarization (*Lp*) factor from each observed peak intensity:

$$Lp = \frac{1 + \cos^2 2\alpha \cdot \cos^2 2\theta}{4\cos\theta \sin^2\theta \cdot (1 + \cos^2 2\alpha)}$$


where α is the diffraction angle of the monochromator.

Note that Equation (23) refers to Bragg–Brentano geometry.


4. Removal of the contribution of the *Lp* factor from the measured intensities *via*:

$$I_{\rm meas}' = \frac{I_{\rm meas}}{Lp}$$


Example: growth of YbPtBi with partial element substitution




Example: growth of YbPtBi with partial element substitution



Only YbPt_{0.9}Au_{0.1}Bi was grown successfully.




Preparation of samples with same stoichiometry is reproducible.



Preparation of samples with varying stoichiometry seems also successful.

N. Ni: PM871-Ex3b



Use of "inner" standard a MUST.

Position of Bragg reflections in powder pattern

$$\lambda = 2d_{hkl}\sin\theta - \frac{1}{d^2} = \frac{1}{V^2} \left(S_{11}h^2 + S_{22}k^2 + S_{33}l^2 + 2S_{12}hk + 2S_{13}hl + 2S_{23}kl \right)$$

$$V = abc\sqrt{1 - \cos^2\alpha - \cos^2\beta - \cos^2\gamma + 2\cos\alpha\cos\beta\cos\gamma}$$

 $S_{11} = b^2 c^2 \sin^2 \alpha$ $S_{22} = a^2 c^2 \sin^2 \beta$ $S_{33} = a^2 b^2 \sin^2 \gamma$ $S_{12} = abc^2 (\cos \alpha \cos \beta - \cos \gamma)$ $S_{13} = ab^2 c(\cos \gamma \cos \alpha - \cos \beta)$ $S_{22} = a^2 b c(\cos \beta \cos \gamma - \cos \alpha)$

Factors affecting peak positions:

$$\Delta 2\theta = \frac{p_1}{\tan 2\theta} + \frac{p_2}{\sin 2\theta} + \frac{p_3}{\tan \theta} + p_4 \sin 2\theta + p_5 \cos \theta + p_6$$


Asymmetry:
$$p_1 = -\frac{h^2 K_1}{3R^2}$$
; $p_2 = -\frac{h^2 K_2}{3R^2}$


In-plane divergence:
$$p_3 = -\frac{\alpha^2}{K_3}$$

Transparency:
$$p_4 = \frac{1}{2\mu_{or}R}$$

Sample displacement:
$$p_5 = -\frac{2s}{R}$$


Zero shift:
$$p_6$$






Combined analysis of series of Bragg reflections (main phase + standard) necessary. $_{\text{N.N.PMST-Exce}}$

Combined fitting of Bragg reflections



Profile of Bragg reflections in powder pattern



Profile of Bragg reflections in powder pattern





Lorentzian $L\!=\!l_0\!\left(1\!+\!\left(\!\frac{2\Theta\!-\!2\Theta_0}{\omega}\!\right)^2\right)^{\!-n}\,n\!=\!1;1.5;2$


Pseudo-Voigt $V = \eta L + (1 - \eta)G$ $(0 \le \eta \le 1)$


Pearson VII $P = l_0 \bigg[1 + \bigg(\frac{2\Theta - 2\Theta_0}{ma^2} \bigg)^2 \bigg]^{-m}$


Parameter: $2\Theta_0$ = peak position; I_0 = peak intensity; ω = FWHM/2; m = shape parameter.











R-values (residuals) - reliability criteria for refinements

$$R = \frac{\sum_{i} |y_{i}(\text{obs}) - y_{i}(\text{calc})|}{\sum_{i} |y_{i}(\text{obs})|}$$

$$R_{wp} = \frac{\sqrt{\sum_{i} w_{i} (y_{i}(\text{obs}) - y_{i}(\text{calc}))^{2}}}{\sqrt{\sum_{i} w_{i} y_{i}(\text{obs})^{2}}} \qquad R_{\exp} = \frac{\sqrt{n - m}}{\sqrt{\sum_{i} w_{i} y_{i}(\text{obs})^{2}}}$$

$$R_{\rm exp} = \frac{\sqrt{n-m}}{\sqrt{\sum_{i} w_i y_i (\rm obs)^2}}$$

$$M = \left(\sum_{i} w_{i} (y_{i}(obs) - y_{i}(calc))^{2}\right)^{1/2}$$

$$R_{Bragg} = \frac{\sum_{j} \left| I_{j}(\text{obs}) - I_{j}(\text{calc}) \right|}{\sum_{j} \left| I_{j}(\text{obs}) \right|}$$




FIG. 2. (Color online) Unit-cell volume and Co concentration determined from WDS measurement as a function of nominal Co concentration. Multiple WDS data points were collected for each nominal x and are each plotted, giving a sense of measured variation in Co concentration. Inset: picture of a representative single crystal over a millimeter grid.






FIG. 3. Unit cell parameters, a and c, as well as unit-cell volume, V, normalized to a_0 =3.9621 Å, c_0 =13.0178 Å, and V_0 =204.3565 ų of undoped BaFe₂As₂ as a function of measured concentration of Co. $x_{\rm WPS}$.

Realized stoichiometry by WDS study; Vegard's law for lattice parameter








High-temperature x-ray diffraction with 2-dimensional detector

Rietveld refinement of powder pattern



Intensity of Bragg reflections in powder pattern

$$I_{hkl} = K \times p_{hkl} \times L_{\theta} \times P_{\theta} \times A_{\theta} \times T_{hkl} \times E_{hkl} \times |F_{hkl}|^{2}$$

Structure amplitude:

$$\mathbf{F}(\mathbf{h}) = \sum_{j=1}^{n} g^{j} t^{j}(s) f^{j}(s) \exp(2\pi \mathbf{i} \mathbf{h} \cdot \mathbf{x}^{j})$$

K = scale factor $p_{hkl} = \text{multiplicity factor}$ $L_{\theta} = \text{Lorentz factor}$

 P_{θ} = polarization factor

 A_{θ} = absorption factor

 T_{hkl} = preferred orientation factor

 E_{hkl} = extinction factor

 g^{j} = population t^{j} = Temperature factor

Lorentz and Polarization factor


$$L_{\theta}$$
 = Lorentz factor

$$L = \frac{1}{\cos\theta \sin^2\theta}$$

 P_{θ} = polarization factor


$$P \propto \frac{1 - K + K \cdot \cos^2 2\theta \cdot \cos^2 2\theta_M}{2}$$

$$LP = \frac{1 + \cos^2 2\theta \cos^2 2\theta_M}{\cos \theta \cdot \sin^2 \theta}, K=0.5$$



Absorption factor

$$A_{\theta}$$
 = absorption factor



$$A = \frac{1}{V} \int_{V} \exp(-\mu_{eff} l) dV$$


$$A = \frac{\mu_{eff}}{2} = \text{const}$$
 (flat opaque sample)

$$A = \frac{1 - \exp(-2\mu_{eff}t/\sin\theta)}{2\mu} \propto 1 - \exp(-2\mu_{eff}t/\sin\theta)$$

flat semitransparent sample


Preferred orientation factor

 T_{hkl} = preferred orientation factor




Extinction factor

 E_{hkl} = extinction factor




Temperature factor


$$t^{j} = \exp\left(-B^{j} \frac{\sin^{2} \theta}{\lambda^{2}}\right)$$
 (isotropic)





Atomic scattering factor

 f^{j} = atomic scattering factor (radial distribution of electrons for x-rays), constant for neutrons











Multi-phase Rietveld refinement





Unit-cell volume of RFeAs(O/F) phase determined by Rietveld analysis

Why and when powder for structure determination?

