## **Neutron and x-ray sources**




Rob McQueeney



Physics 590

## Different probes










	NEUTRONS	X-RAYS	ELECTRONS
Wavelength range	0.4 - 10 Å	0.1 - 5 Å	0 .04 - 0.2 Å
Energy range	0.001 - 0.5 eV	3000 - 100000 eV	6000 - 120000 eV
Cross-section	10 <sup>-25</sup> barns	10 <sup>-25</sup> Z <sup>2</sup> barns	~10 <sup>-22</sup> barns
Penetration depth	~ cm	~ μ <b>m</b>	~ nm
Typical flux	10 <sup>11</sup> s <sup>-1</sup> m <sup>-2</sup>	10 <sup>24</sup> s <sup>-1</sup> m <sup>-2</sup>	10 <sup>26</sup> s <sup>-1</sup> m <sup>-2</sup>
Beam size	mm-cm	μ <b>m-mm</b>	nm-μm
Typical sample	Any bulk sample	Small crystals, powders, surfaces	Surfaces, thin films, grains, gases
Techniques	Diffraction Inelastic scattering Reflectivity	Diffraction Photon absorption Photoemission Inelastic scattering	Microscopy Diffraction Emission spectroscopy EELS
Phenomena	Magnetic/crystal structures collective excitations (phonons, spin waves) electronic excitations (crystal- field, spin-orbit)	Crystal structures, electronic transitions (photoemission, absorption),	microstructure crystal structures electronic transitions

Physics 590

## **Cross-sections**





### Neutrons


- Random with Z
- Depends on isotope
- Depends on nuclear spin
- Absorption can be problem

		Abundance	Cross-	Absorption
		(%)	section (bn)	(bn)
	Gd		180	49700
	152Gd	0.2	13	735
	154Gd	2.1	13	85
	155Gd	14.8	66	61100
	156Gd	20.6	5	1.5
	157Gd	15.7	1044	259000
	158Gd	24.8	10	2.2
Phy	160Gd	21.8	10.52	0.77
	·	•	·	

## **Producing neutrons**




**Fission**Nuclear reactor



Moderators →Cold-Thermal

# **Spallation**Particle accelerator



Moderators →Cold-Epithermal


## **Neutrons by reactor fission**






High flux isotope reactor - ORNL








Physics 590








### High Flux Isotope Reactor at Oak Ridge National Laboratory



ORY

#### The United States' highest flux reactor-based source of neutrons for condensed matter research



07-G00244E/arm

Under consideration





## Neutrons by pulsed spallation







### **Spallation Neutron Source (ORNL)**

## **Target-moderator system**





SNS liquid Hg target





Fig. 3. Horizontal cross-section of the flux-trap moderators.



### Spallation Neutron Source at Oak Ridge National Laboratory

#### The world's most intense pulsed, accelerator-based neutron source

#### Backscattering Spectrometer (BASIS) . BL-2

Dynamics of macromolecules, constrained molecular systems, polymers, biology chemistry, materials science

**Spallation Neutrons and Pressure** 

Diffractometer (SNAP) • BL-3

Materials science, geology, earth and

environmental sciences

Chris Tulk + 865.576.7028 + tulkca@ornl.go

Eugene Mamontov · 865.574-5109 · mamontove@ornl.gov

#### Nanoscale-Ordered Materials Diffractometer (NOMAD) • BL-1B (2010)

Liquids, solutions, glasses, polymers, nanocrystalline and partially ordered complex materials

Joerg Neuefeind • 865.241.1635 • neuefeindic@ornl.gov

#### **Wide Angular-Range** Chopper Spectrometer (ARCS) • BL - 18

Atomic-level dynamics in materials science, chemistry, condensed matter sciences

Doug Abernathy · 865.576.5105 · abernathydl@ornl.gov

#### Fine-Resolution Fermi Chopper Spectrometer (SEQUOIA) • BL - 17 (2008)

Dynamics of complex fluids, quantum fluids, magnetism, condensed matter, materials science

Garrett Granroth • 865.576.0900 • granrothge@ornl.gov

#### **Ultra-Small-Angle Neutron Scattering Instrument** (TOF-USANS) • BL-1A (2012\*

Life sciences, polymers, materials science. earth and environmental sciences Michael Agamalian · 865.576.0903 ·

#### Chemical Spectrometer (VISION) • BL-16B (2011)

Vibrational dynamics in molecular systems, chemistry Christoph Wildgruber • 865.574.5378 • wildgrubercu@ornl.gov

BL - 16A

#### **Neutron Spin Echo Spectrometer** (NSE) • BL - 15 (2009)

High-resolution dynamics of slow processes, polymers, biological macromolecules Michael Ohl • 865.574.8426 • ohlme@ornl.gov

#### **Hybrid Spectrometer** (HYSPEC) • BL-14B (2011)

Atomic-level dynamics in single crystals, magnetism, condensed matter sciences

Mark Hagen • 865.241.9782 •

BL-14A

#### Magnetism Reflectometer • BL-4A

Chemistry, magnetism of layered systems and interfaces Valeria Lauter • 865.576.5389 • lauterv@ornl.gov

#### Liquids Reflectometer • BL-4B

Interfaces in complex fluids, polymers, chemistry John Ankner • 865.576.5122 • anknerjf@ornl.gov

#### **Cold Neutron Chopper** Spectrometer (CNCS) • BL - 5

Condensed matter physics, materials science, chemistry, biology, environmental science

Georg Ehlers • 865.576.3511 • ehlersg@ornl.gov

#### **Extended Q-Range Small-Angle Neutron** Scattering Diffractometer (EQ-SANS) • BL-6 (2008)

Life science, polymer and colloidal systems, materials science. earth and environmental sciences

Jinkui Zhao · 864.574.0411 · zhaoj@ornl.gov

**Elastic Diffuse Scattering** Spectrometer (CORELLI) . BL-9 (2013)

> Detailed studies of disorder in crystalline materials Feng Ye • 865.576.0931 • yef1@ornl.gov

**BL-10** 

#### Macromolecular Neutron Diffractometer (MaNDi) • BL-11B (2012)

Atomic-level structures of membrane proteins, drug complexes, DNA

Leighton Coates

#### **Fundamental Neutron** Physics Beam Line • BL-13 (2008)

Fundamental properties of neutrons

Geoffrey Greene • 865.574.8435 • greenegl@ornl.gov

#### **Single-Crystal Diffractometer** [TOPAZ] • BL-12 (2009)

Atomic-level structures in chemistry, biology, earth science, materials science, condensed matter physics

> Christina Hoffmann - 865.576.5127 hoffmanncm@ornl.gov

Powder Diffractometer (POWGEN) . BL-11A (2008)

Atomic-level structures in magnetism, chemistry, materials sciences Jason Hodges • 865.576.7034 • hodgesj@ornl.gov

#### \* Scheduled commissioning date



**Engineering Materials Diffractometer** (VULCAN) • BL-7 (2008)

BL-8A

BL-8B

Mechanical behaviors, materials science, materials processing

Xun-Li Wang • 865.574.9164 • wangxl@ornl.gov










## **Powder diffraction**




### **Determine the crystal structure**



## **TOF** powder diffraction







POWGEN @ SNS

Time-of-flight  $\tau = L/v = \lambda mL/h = 2mLd\sin\theta/h$ 

Physics 590






Time, wavelength, or d-spacing






# Single-crystal: more detail than powders Wide angle diffraction: Get an overview of everything















## Single crystal diffraction




### Triple-axis diffraction: focus in on specific points of interest



HB-1A 3-axis spectrometer






Orbital ordering in YVO<sub>3</sub>


## **Production of x-rays**





X-ray tube

Rotating anode






## How do we get x-rays?



- Static charge ---- Electric field
- Charge moving at constant v ---- Magnetic field

Accelerating charge --- Electromagnetic radiation



Synchrotron radiation is

highly collimated

highly linearly polarized

highly brilliant


continuous wavelength distribution (beyond Cu, Mo, etc..)



## **Advanced Photon Source**

### **Synchrotron**








## **Insertion Devices**











## Main Undulator Line





## Places to go





### **Further references**



### General neutron scattering

- G. Squires, "Intro to theory of thermal neutron scattering", Dover, 1978.
- S. Lovesey, "Theory of neutron scattering from condensed matter", Oxford, 1984.
- R. Pynn, <a href="http://www.mrl.ucsb.edu/~pynn/">http://www.mrl.ucsb.edu/~pynn/</a>.

### Structural refinements

- GSAS <a href="http://www.ncnr.nist.gov/xtal/software/gsas.html">http://www.ncnr.nist.gov/xtal/software/gsas.html</a>
- FullProf http://www.ill.eu/sites/fullprof/

### How to get beam time

- Talk to one of us at Ames about your experiment
- We can identify a suitable instrument
- Talk to instrument scientist
- Write a beamtime request